
# Boom Stand B Auslegerstativ B

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English

Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

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Issued by: Carl Zeiss Microscopy GmbH  
Carl-Zeiss-Promenade 10  
07745 Jena, Germany  
  
microscopy@zeiss.com  
www.zeiss.com/microscopy

Carl Zeiss Microscopy GmbH  
Königsallee 9-21  
37081 Göttingen, Germany

**Number of this manual:** 435424-7144-010  
**Date of issue:** Version 3, 12/19/2016

**CONTENTS**

	Page
<b>1 INTRODUCTION .....</b>	<b>2</b>
1.1 General information.....	2
1.2 Notes on instrument safety .....	2
1.3 Notes on warranty .....	3
<b>2 DESCRIPTION .....</b>	<b>4</b>
2.1 Intended use.....	4
2.2 System overview .....	4
2.3 Components and controls .....	5
2.4 Technical data .....	6
2.5 Dimensions and workspaces.....	7
<b>3 UNPACKING AND ASSEMBLY .....</b>	<b>9</b>
3.1 Unpacking and assembling the Boom Stand B .....	9
3.2 Installing Stemi mount and microscope body .....	11
3.3 Mounting and adjusting external illumination on the Stemi .....	11
<b>4 OPERATION .....</b>	<b>12</b>
4.1 Vertical adjustment on the vertical column .....	12
4.2 Horizontal adjustment of the boom arm at the cross piece .....	13
4.3 Boom rotation about the vertical column.....	13
4.4 Tilting on the Tilting Head .....	14
4.5 Rotating the Stemi mount .....	14
4.6 Focusing drive on the Stemi mount .....	14
4.7 Application examples .....	15
4.7.1 Vertical viewing position .....	15
4.7.2 Oblique viewing position.....	15
4.7.3 Observation position for maximum sample dimensions .....	16
4.7.4 Horizontal viewing position .....	16
<b>5 DISPOSAL OF THE PRODUCT .....</b>	<b>17</b>

## 1 INTRODUCTION

### 1.1 General information

The Boom Stand B is an accessory for the Stemi 305 or Stemi 508 stereo microscopes. It can also be used with SteREO Discovery configurations whose weight does not exceed the working load of the Stemi mount (5 kg).

This manual covers installation and operation of the Boom Stand B.

For first-time set-up of the stereo microscopes, also observe the separate operating manuals for Stemi 305, Stemi 508 or Stemi SteREO Discovery and the operating manuals of the external cold light sources.

The following warning and information symbols are used in this operating manual:

#### CAUTION

*This symbol indicates a potential hazard to the user.*

#### ATTENTION

*This symbol indicates a potential hazard to the instrument or system.*

#### ATTENTION

*Read the operating manual.*

#### NOTE

*This symbol designates information that should be closely followed.*

### 1.2 Notes on instrument safety

The Boom Stand B for stereo microscopes, including original accessories, may only be used for the applications described in this operating manual. The manufacturer cannot assume any liability for other applications of the instrument, including those of individual modules or single parts. This also applies to any service or repair work that is not carried out by authorized service personnel. Failure to comply with this shall render all warranty claims invalid.

The safety instructions in the operating manuals of all components used with the Boom Stand B, such as the microscope body, illumination equipment etc., must be observed.

---

### 1.3 Notes on warranty

The manufacturer guarantees that the instrument is free of material or manufacturing defects when delivered. Possible defects must be notified to us immediately and steps taken to minimize damage. If notified of such a defect, the manufacturer is obligated to rectify it at his discretion, either by repairing the instrument or delivering an intact replacement. No guarantee is provided for defects caused by natural wear (wearing parts in particular) and improper use.

The instrument manufacturer shall not be liable for damage caused by faulty operation, negligence or any other tampering with the instrument, particularly the removal or replacement of instrument components, or the use of accessories from other manufacturers. Such actions will render any warranty claims invalid.


With the exception of the work described in this manual, no maintenance or repair work is to be carried out on these microscopes. Repairs may only be performed by ZEISS Service or individuals specially authorized by ZEISS Service. In the event of a problem with the instrument, please contact the Carl Zeiss microscopy service team or your local ZEISS overseas representative.

**2 DESCRIPTION****2.1 Intended use**


The Boom Stand B is an accessory for the Stemi 305 or Stemi 508 stereo microscopes and SteREO Discovery.

Due to its long boom arm and 650 mm high stand column, the boom stand allows the investigation of particularly high or large-scale samples, or the observation of many small samples in a large sample compartment.

The Boom Stand B is intended for use in biological and veterinary laboratories as well as in industrial production and quality assurance.

**2.2 System overview**

2.3 Components and controls

- 1 Star knob for horizontal adjustment
- 2 Cross piece
- 3 Star knob for vertical adjustment
- 4 Wing screw for safety ring
- 5 Safety ring
- 6 Adapter BMS
- 7 Clamping lever for tilting
- 8 Tilting Head
- 9 Star knob on Stemi mount
- 10 Focus drive on Stemi mount
- 11 Clamping screw on Stemi mount for microscope body

**Fig. 1 Components and controls on the Boom Stand B or Stemi mount**

English

**2.4 Technical data****Dimensions and weight:**

Boom Stand B, complete equipment with Stemi mount .....	approx. 25.3 kg
consisting of:	
- Stand Base B 32 (280 mm x 200 mm x 47 mm) .....	approx. 20.5 kg
- Column 32/650 mm (Ø 32 x 650 mm) .....	approx. 1.4 kg
- Column 32/450 mm (Ø 32 x 450 mm) .....	approx. 1.1 kg
- Cross piece B 32 .....	approx. 0.5 kg
- Adapter B-32 BMS 90° .....	approx. 0.2 kg
- Tilting Head B 32 .....	approx. 0.7 kg
- Safety ring B 32 .....	approx. 0.1 kg
- Stemi mount .....	approx. 0.8 kg

**Adjustment ranges:**

horizontal on cross piece B 32 .....	max. 380 mm
vertical on column 32/650 mm (without safety ring) .....	max. 560 mm
rotatable on column 32/650 mm .....	360°

**Accessible object positions:**

Swivel radius about vertical column, measured up to

axis of the Tilting Head .....	75 to 455 mm
optical axis of the microscope body with radially oriented Stemi mount .....	230 to 610 mm

2.5 Dimensions and workspaces

English


Fig. 2 Dimensions and workspaces of Boom Stand B with microscope body


Fig. 3 Dimensions and working positions of Stand Base B 32

Item No.	Dimensions / Workspace	Remarks
1	up to 380 mm	Horizontal adjustment of boom arm
2	up to 560 mm	Vertical adjustment of stand column
3	up to 455 mm	Distance between stand column and Tilting Head axis
4	up to 610 mm	Distance between stand column and optical axis of the microscope body
5	up to 555 mm	Height of Stemi mount (microscope support) above work surface – Tilting Head mounted downwards
6	up to 865 mm	Height of Stemi mount (microscope support) above work surface – Tilting Head mounted upwards
7a	650 mm	Minimum depth of the work surface
7b	675 mm	Minimum depth of the work surface
7c	850 mm	Minimum depth of the work surface
8a	up to 370 mm	Distance between center of the object field and Stand Base
8b	up to 475 mm	Distance between center of the object field and Stand Base
8c	up to 570 mm	Distance between center of the object field and Stand Base

### 3 UNPACKING AND ASSEMBLY

#### 3.1 Unpacking and assembling the Boom Stand B

The Boom Stand B consists of the following components:


- Stand Base B 32
- Column 32/650 (vertical)
- Column 32/450 (horizontal)
- Safety ring B 32
- Cross piece B 32
- Adapter B 32-BMS 90°
- Tilting Head B 32/115 mm with end piece


**CAUTION**

When unpacking and installing the Stand Base B 32, consider that it has a weight of approx. 20 kg.


- Remove all components from the packaging (verify that all components are present).
- Put the Stand Base (Fig. 4/1) in the desired position on a stable surface.
- Remove the two clamping screws 3 mm (Fig. 4/2) on the Stand Base using the enclosed Allen screwdriver.
- Fit the vertical column (Fig. 4/3) as far as it will go into the Stand Base and clamp using the two clamping screws.
- Slip the safety ring (Fig. 4/4) over the column and clamp it at the desired height using the wing screw (Fig. 4/5).
- Unscrew both star knobs (Fig. 5/2 and 5) slightly from the cross piece until they do not project into the respective hole.
- Slip the cross piece (Fig. 5/1) over the vertical column (Fig. 5/3) up to the protrusion of the safety ring and clamp in the desired position using the star knob (Fig. 5/2).
- Push the horizontal column (Fig. 5/4) into the hole of the cross piece and clamp it in the desired horizontal position using the star knob (Fig. 5/5).

**Fig. 4** Mounting the vertical column 32/650 onto the Stand Base B 32


**Fig. 5** Mounting the horizontal column 32/450

**Fig. 6 Mounting the adapter BMS and Tilting Head**

- Remove the two clamping screws 3 mm (Fig. 6/3) on the adapter BMS (Fig. 6/2) using the enclosed Allen screwdriver.
- Fit the adapter BMS to the end of the horizontal column (Fig. 6/1) and tighten both clamping screws 3 mm.

**Fig. 7 Mounting the Tilting Head**

- Loosen the two clamping screws 3 mm (Fig. 7/3) in adapter BMS and pull the plastic clamp forward.
- Fit the Tilting Head (Fig. 7/2) into the receptacle of the adapter BMS (Fig. 7/1), press the plastic clamp against the Tilting Head and tighten both clamping screws (right and left equally).
- If necessary, align the Tilting Head perpendicularly by rotating the horizontal column. To do this, first loosen the star knob (Fig. 5/5) and then tighten again.

Depending on the weight of the components to be attached to Boom Stand B (e.g. Stemi mount, microscope, illuminator), position and clamp the boom vertically and horizontally for maximum space for the sample. Always check and ensure the stability.

After assembly of all components, check the stability and reduce the boom length or height as necessary, or change the position of the boom relative to the Stand Base (see also Fig. 3 page 8).

### 3.2 Installing Stemi mount and microscope body

- Loosen the clamping lever (Fig. 8/1) on the Tilting Head (Fig. 8/2).
- Align Tilting Head (Fig. 8/2) vertically or set in desired angular position, then re-tighten the clamping lever.

 The angular position of the clamping lever (Fig. 8/1) can be set to a position that is favorable for the clamping procedure. To this end, pull the lever until it comes out of the click stop and turn it into the desired angular position. After releasing the lever, it will automatically snap in again by a spring mechanism.

- Unscrew the end piece (Fig. 8/5) from the Tilting Head (Fig. 8/2).
- Slip Stemi mount (Fig. 8/3) from the bottom over the Tilting Head and secure it in the desired position by means of the star knob (Fig. 8/4).
- Screw the end piece back into Tilting Head.

#### ATTENTION

Before inserting the microscope body, make sure that its objective does not collide with the Stand Base or the mounting surface of the stand when inserting it. If necessary, change the vertical adjustment.

- Insert the microscope body, e.g. Stemi 508 (Fig. 8/7), into the mounting hole of the Stemi mount and fix using the clamping screw (Fig. 8/6). Check the stability of the stand and reduce the boom length or height, if necessary, or change the position of the boom to the Stand Base (see Fig. 3).
- If present, screw a front optics or objective.

 If the focusing drive of the Stemi mount drops unintentionally, the torque of the focusing drive must be adjusted to the weight of the microscope. To do this, rotate the left and right focus drive in opposite directions.


Fig. 8 Installing Stemi mount and microscope body

### 3.3 Mounting and adjusting external illumination on the Stemi

The external illumination (fiber-optic or LED) should be mounted and operated on the microscope body used according to the respective operating manual.

## 4 OPERATION

The Boom Stand B allows to adjust the stereo microscope optimally to large-size or voluminous objects. The following adjustments are possible:

- Vertical adjustment of the vertical column, see Section 4.1
- Horizontal adjustment of the horizontal column, see Section 4.2
- Boom rotation about the vertical column, see Section 4.3
- Tilting at the Tilting Head, see Section 4.4
- Rotation of the Stemi mount, see Section 4.5

### CAUTION

The safety ring (Fig. 1/5) must always be clamped directly beneath the cross piece (Fig. 1/2)! The clamps of safety ring (Fig. 1/4) and cross piece (Fig. 1/3) may only be released alternately; otherwise, the cross piece may inadvertently slip down together with the complete structure, risking injury to the operator or damage to the microscope and the object to be examined.

### 4.1 Vertical adjustment on the vertical column

#### Vertical adjustment, upward:

- Make sure that the safety ring (Fig. 1/5) is clamped. If not, tighten the wing screw (Fig. 1/4).
- Loosen star knob (Fig. 1/3) and push the cross piece (Fig. 1/2) with the complete structure upwards into the desired position.
- Retighten the star knob (Fig. 1/3).
- Hold the safety ring (Fig. 1/5) tight, release wing screw (Fig. 1/4) and push the safety ring up until it touches the cross piece.
- Retighten the wing screw (Fig. 1/4).

#### Vertical adjustment, downward:

### ATTENTION

Before any downward vertical adjustment, check how far the boom arm can be moved downward without the microscope colliding with the mounting surface, base plate or object to be examined.

- Make sure the cross piece is clamped; if not, tighten the star knob (Fig. 1/3).
- Release the wing screw (Fig. 1/4) and push safety ring (Fig. 1/5) down into the desired position.
- Retighten the wing screw (Fig. 1/4).

- Hold the cross piece (Fig. 1/2) with complete structure and release the star knob (Fig. 1/3).
- Slowly lower the cross piece (Fig. 1/2) until it rests on the safety ring.
- Retighten the star knob (Fig. 1/3).

#### 4.2 Horizontal adjustment of the boom arm at the cross piece

##### CAUTION

The stability of the stand is no longer ensured while the horizontal position of the boom arm is being changed.

The horizontal position of the boom arm, alignment of the Stand Base and weight of the microscope structure must be coordinated to ensure stability. Otherwise, the stand may tip over, with a risk of injury to the operator or damage the microscope or object to be examined.

- Loosen star knob (Fig. 1/1).
- Extend or retract the horizontal column with the complete structure to the desired length.
- Retighten the star knob (Fig. 1/1).

#### 4.3 Boom rotation about the vertical column

##### CAUTION

The stability of the stand is no longer ensured while the horizontal position of the boom arm is being changed. It is therefore important to rotate only slowly and watch for stability. If necessary, reduce the horizontal position of the boom arm.

Take care to prevent any collision with persons, the object to be examined or any other surrounding objects. Take into account the length and routing of the optical fiber.

The safety ring (Fig. 1/5) must always be clamped directly beneath the cross piece (Fig. 1/2)! The clamps of safety ring (Fig. 1/4) and cross piece (Fig. 1/3) may only be released alternately; otherwise, the cross piece may inadvertently slip down together with the complete structure, risking injury to the operator or damage to the microscope and the object to be examined.

- Push the safety ring (Fig. 1/5) under the cross piece (Fig. 1/2) and clamp it tightly (Fig. 1/4).
- Loosen the star knob Fig. 1/3) to rotate about the vertical column and swivel the cross piece with the complete structure.
- Retighten the star knob (Fig. 1/3).

#### 4.4 Tilting on the Tilting Head

##### ATTENTION

Before tilting, check how far the Tilting Head can be inclined without the microscope colliding with the double boom, the object to be examined or the mounting surface or Stand Base. If necessary, increase the height of the boom arm.

- Loosen the clamping lever (Fig. 1/7) on the Tilting Head Fig. 1/8).
- Incline the Tilting Head (Fig. 1/8) with Stemi mount and microscope into the desired position.
- Retighten clamp lever (Fig. 1/7).

The microscope can also be tilted, if necessary, by turning the horizontal column in the cross piece, for example, when setting the horizontal viewing position, see Section 4.7.4. To do this, loosen the star knob (Fig. 1/1) and then tighten again.

#### 4.5 Rotating the Stemi mount

##### ATTENTION

Before rotating the Stemi mount, make sure that the end piece on the Tilting Head is screwed on and check how far it can be rotated without the microscope colliding with the structure.

- Loosen star knob (Fig. 1/9) on Stemi mount.
- Rotate Stemi mount and microscope into the desired position.
- Retighten the star knob (Fig. 1/9).

If necessary, the microscope body also be rotated within the Stemi mount, e.g. when setting a horizontal viewing position, see Section 4.7.4. Prior to rotating the microscope body, loosen the clamping screw (Fig. 1/11) and then retighten.

#### 4.6 Focusing drive on the Stemi mount

- Turn the focusing drive on the Stemi mount (Fig. 1/10) to adjust the height of the microscope or to focus on the object.
- By reverse rotation of the left and right focusing drive, the torque of the focusing drive can be adapted to the needs of the user.


##### ATTENTION

The movement should not be set too freely to prevent the drive from inadvertent slipping. This could lead to damage of the microscope or the object.


**4.7 Application examples**

**4.7.1 Vertical viewing position**


**4.7.2 Oblique viewing position**

Enables the variation of viewing angle and viewing height for an ergonomic seat position and relaxed observation.

### 4.7.3 Observation position for maximum sample dimensions


#### CAUTION

In this position, take care that the clamping lever (Fig. 1/7) of the Tilting Head and star knob (Fig. 1/1) are firmly tightened. If they loosen inadvertently, there is the risk of injury to the operator or damage to the microscope.

In this observation position the viewing angle and height can also be varied by tilting the Tilting Head.

### 4.7.4 Horizontal viewing position

When observing upright samples such as sculptures or paintings in the field of restoration.

---

## 5 DISPOSAL OF THE PRODUCT

The product was developed, tested and manufactured in conformance with the applicable regulations and environmental directives of the European Union.

The product and accessories satisfy the requirements of EU Directives 2011/65/EC (RoHS) and 2012/19/EC (WEEE), as well as the German Act on Electronic and Electrical Equipment (ElektroG).

For details of disposal and recycling, please contact your sales partner or ZEISS Service.

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