
Stand M LED
Stativ M LED
Statif MLED
Estativo M LED
Tripé M LED
Stativo M LED
Штатив M LED
M LED 鏡臂
スタンド M LED
M LED الحامل الثلاثي

Operating manual
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo
Manual de instruções
Istruzioni per l'uso
Инструкция по применению
操作手册
取扱説明書
دليل المستخدم

ZEISS

English

Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

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Number of this manual: 435425-7144-008
Date of issue: Version 1, 06/30/2016

1 INTRODUCTION

Stand M LED is an accessory for the Stemi 305 or Stemi 508 stereo microscopes.

The present operating manual is valid only in conjunction with the operating manuals of Stemi 305 or Stemi 508. Users must read them carefully, as well as the instructions regarding their safety and installation aspects. The safety information concerning the Stand K EDU/LAB that is contained, also applies for Stand M LED.

Please find actual and complete operation manuals in several languages in the Download Area of the following Websites:

www.zeiss.de/Stemi305

www.zeiss.com/Stemi305

www.zeiss.de/Stemi508

www.zeiss.com/Stemi508

Stand M LED is the large desktop stand for the Stemi 305 or Stemi 508 stereo microscopes. Thanks to its generously proportioned work surface and 350 mm high focusing column it is suitable both for use with elongated or tall samples as well as for parallel working with multiple samples.


Stand M LED contains electronics for controlling the K/M LED reflected and transmitted light illumination units (accessories).

This operating manual describes the installation and operation of Stand M LED, as well as the optional transmitted light illuminators:

- Transilluminator Brightfield M LED
- Transillumination Unit M LED.

For detailed information on reflected-light illuminators refer to the Stemi 305 / 508 operating manuals.

2 MICROSCOPE SYSTEM

- | | |
|---|--|
| <p>1 Eyepiece in eyepiece tube (e.g. eyepiece 10x/23)</p> <p>2 Microscope body Stemi 508</p> <p>3 Zoom button for magnification settings, including switchable click stops</p> <p>4 Stand M LED</p> <p>5 Focus drive for focusing the sample</p> <p>6 Memory keys for illumination settings</p> <p>7 Rotary pushbutton to control the reflected light illuminators K LED (On/Off/Dimming)</p> | <p>8 Rotary pushbutton to control the transmitted light illuminators M LED (On/Off/Dimming)</p> <p>9 Reflected light illuminator (e.g. K LED spot light)</p> <p>10 Insert plate for sample positioning</p> <p>11 Lever to switch transmitted light illumination - brightfield/darkfield - on Transilluminator Brightfield</p> <p>12 Rotary / slider button to set transmitted light illumination - brightfield, oblique light or darkfield - on Transillumination Unit M LED</p> |
|---|--|

Fig 1 Stand M LED with Transilluminator Brightfield M LED (left) and Transillumination Unit M LED

3 MOUNTING THE TRANSMITTED LIGHT ILLUMINATORS M LED

CAUTION

Switch off the stand and unplug the power supply cable from the stand before mounting the transmitted light illuminator.

3.1 Mounting the Transilluminator Brightfield M LED

- Unplug the power supply cable from the stand.
- Gently tilt the stand backwards and put it down.
- Unscrew 3 mm Allen screw (Fig. 2/1) and put it to one side.
- Remove cover flaps (Fig. 2/2 and 3).
- Insert the Transilluminator Brightfield so that the two pins (Fig. 2/4) click into the recesses (Fig. 2/5) on the bottom of the stand.
- Press the Transilluminator Brightfield against the stand base and hold tight.
- Screw the 3 mm Allen screw into the tapped hole (Fig. 2/6) of the stand until the tip of the screw engages in the hole of the Transilluminator Brightfield.
- Tighten the Allen screw hand-tight.
- Return the stand to upright and, if necessary, prevent wobbling by adjusting the height-adjustable foot of the stand base (Fig. 2/7).
- Plug the power supply cable back into the stand.


Fig. 2 Mounting the Transilluminator Brightfield M LED


Fig. 3 Mounting the Transillumination Unit

3.2 Mounting the Transillumination Unit M LED

- Unplug the power supply cable from the stand.
- Place the Transillumination Unit on a stable surface.
- Unscrew the 3 mm Allen screw (Fig. 3/1) on the front of the stand and place it to one side.
- Screw the height-adjustable stand foot (Fig. 3/2) in completely.
- Remove the cover flap (Fig. 3/3).
- Carefully mount the stand onto the Transillumination Unit from the rear. Push it forwards until the stand plate and transmitted light housing are flush. The two lugs (Fig. 3/5) must slot into the recesses (Fig. 3/4).
- Screw the 3 mm Allen screw into the tapped hole (Fig. 3/6) of the Transillumination Unit and tighten by hand.
- Return the stand to upright and, if necessary, prevent wobbling by adjusting the height-adjustable foot of the Transillumination Unit (Fig. 3/7).
- Plug the power supply cable back into the stand.

To disassemble the Transillumination Unit M LED, unscrew the Allen screw, carefully push the stand backwards and remove it.

4 OPERATION

4.1 Operation of Stand M LED with Transilluminator Brightfield

- Press the lower knob (Fig. 4/2) to switch on the transmitted light illuminator and adjust it by turning the illumination intensity control.
- Toggle between brightfield (BF) and all-sided darkfield (DF) illumination by moving the lever (Fig. 4/3).

 Press/turn the upper illumination button (Fig. 4/1) to add a variable amount of reflected light to the set transmitted light.

Detailed information on operating the optional reflected-light illuminators K LED can be found in the Stemi 305 / 508 operating manual.


Fig. 4 Operating the Transilluminator Brightfield

4.2 Operation of Stand M LED with Transillumination Unit

- Press the lower knob (Fig. 5/2) to switch on the transmitted light illuminator and adjust it by turning the illumination intensity control.
- To adjust the transmitted light brightfield illumination, slide the rotary / slider button (Fig. 5/3) to the rear position (BF).
- Then adjust the tilting angle of the rotating mirror by turning the button until the specimen is brightly illuminated from below.

 The rotating mirror has plain and frosted mirror sides: use the plain side for high-contrast brightfield, and the diffuse (frosted) mirror for homogenous brightfield.

- To set lateral transmitted light darkfield illumination, slide the control knob (Fig. 5/3) to the front position (DF).


Fig. 5 Operating the Transillumination Unit

- Then adjust the tilting angle of the rotating mirror by turning the knob until the sample structures appear brightly illuminated against a dark background.

Slide the rotary / slider button (Fig. 5/3) to a central position (RC) and then tilt the mirror to create oblique illumination which reveals the structures of unstained specimens in relief contrast.

Press/turn the upper illumination button (Fig. 5/1) to add a variable amount of reflected light to the set transmitted light.

Detailed information on operating optional reflected-light illuminators K LED can be found in the Stemi 305 / 508 operating manual.


Fig. 6 Operating the memory keys

4.3 Operation of the memory keys

Stand M LED has three memory buttons (Fig. 6/1) to save different illumination settings and to permit rapid switching between different reflected, transmitted and mixed illumination types.

To save an illumination setting, proceed as follows:

- Turn on the desired illumination types (reflected and / or transmitted light) by pressing the top (Fig. 6/2) or bottom knob (Fig. 6/3).
- Adjust the desired intensity of the reflected and / or transmitted light illumination by turning the top (Fig. 6/2) or bottom knob (Fig. 6/3).
- Choose one of the three memory buttons (Fig. 6/1) and hold it down for 3 seconds. All 3 buttons will light up when the settings have been successfully saved.

- Repeat the process to assign the other two memory keys.
- To recall an illumination setting, press the relevant memory key shortly.

The switching status (on / off) and the brightness of all adapted illuminators are saved. Exception: In the case of the K LED ring illumination, only the switching status (on / off) is saved. The brightness and segment settings of the ring illumination are not saved and remain unchanged when memory keys M1-M3 are pressed.

Factory reset: With Stand M LED switched off, press all three memory buttons M1-M3 simultaneously while switching on the stand. Memory slots M1 / 2 / 3 will be emptied and all connected illuminators switched on at 50% brightness. Briefly pressing buttons M1-M3 now causes no change in the illumination until a new setting is made.

5 TECHNICAL DATA

Weight4.6 kg

Optical risk group classification acc. to DIN EN 62471:2009


SystemLED risk group 2 acc. to
DIN EN 62471:2009

Transillumination Unit M LED in Stand M LED LED risk group 2 acc. to
DIN EN 62471:2009


Transilluminator Brightfield M LED in Stand M LED..... LED risk group 2 acc. to
DIN EN 62471:2009

 For detailed information on technical data refer to the operating manual for Stemi 305 / 508.

Dimensions of Stand M LED with Transilluminator Brightfield M LED

Dimensions of Stand M LED with Transillumination Unit M LED

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